



# BLUEBERRIES

Cardboard trays with PE flow packs

**PHOTO CONDITIONS:** *Pictures do not show flow packs to improve the quality of the photo.*

## RESULTS:

**Appearance:** After 20 days at 7°C, blueberries contained in the package treated with cronogard® are more hydrated and do not show mold and rot while blueberries contained in the reference are moldy and dehydrated.

**Consistency:** Blueberries packed in treated trays are still firm, while those in reference packaging are much softer and overripe with evidence of molds.

## TEST CONDITIONS:

20 days at 7°C

## RAW MATERIAL and APPLICATION TECHNIQUE:

cartonboard 100% virgin fibres  
370 g/m<sup>2</sup>; flexo coating 2 layers

## TYPE OF PACKAGING: :

cronogard® treated cardboard trays with untreated PE flow packs

## SAMPLE REFERENCE:

untreated trays and PE flow packs

## FILLER and LACQUER:

CC19-A6B607-CIF

WITH  
CRONOGARD®



WITHOUT  
CRONOGARD®

